

ora-laser 2100

Our therapy laser
for an extra gentle application.

- **No protective goggles necessary!**
- **intuitive operation by integrated Touch-Display**
- **200 preinstalled programs and 100 further program places for your individual programs**
- **exact doses by processor controlled laser emission and applicator identification**
- **up to 100 Milliwatt power**
- **application programs:**
 - o dental
 - o human
 - o cosmetic
 - o acupuncture
- **3 applicators (autoclavable)**
- **1 special applicator (autoclavable) with acupuncture point detection system**



Technical data

Diode laser with 810nm technology, laser class 3R, standard applicators in 3 diameters and 1 special applicator with integrated acupuncture point detection system - all autoclavable, broadband frequency up to 10.000 Hz., handpiece with applicator identification for correct doses, individual pin code entry, multiple Touch-Panel.

ora-laser 2100

applications: medical	applications: dental	applications: cosmetic
Burns 2nd + 3rd degree	Abscesses	Basic Cosmetic Treatment 1
Bursitis	Acidity Trauma	Basic Cosmetic Treatment 2
Carpal-Cannel-Synd.	Allergy to prosthetic metals	Acne 1
Claudicatio Intermittent	Aphthae	Acne 2
Colitis Ulcerous	Bio-stimulation of wounds and cicatrices	Acne 3
Contorsion	Cicatrix Smoothing	Blackheads, Comedos
Contusion	Decubitus	Couperosis (small)
Coxarthrosis	Dentitio Difficilis (difficult dentition)	Couperosis (big)
Decubitus	Dolor Post 1 (post extractional pain)	Irritations (general; small areas of skin)
Diarrhea chronic, non infectious	Dolor Post 2 (post-e. pain, prophylactic)	Irritations (general; large areas)
Distorsion	Fresh Wounds	Irritations (hair removal; large areas)
Dysmenorrhea	Gingivitis	Irritations (hair removal; small areas)
Eczema	Granulomas	Strengthening of connective tissue (on belly, upper arms or thighs)
Epicondylitis Humeri	Hematomas 1	Breast Lifting
Erysipelas	Hematomas 2	Breast Lifting (punctually)
Facialis-parese	Herpes Labialis	Cellulite
Fissures	Jaw Ostitis after Surgery	Facelifting
Fractures	Jaw Trismus	Smoothing of Cicatrices (small, fresh)
Furuncle	Nausea 1	Bigger Cicatrices
Gout	Nausea 2	Cicatrices, old and small
Headache	Neuralgic Disorders 1	Cicatrices, old and big
Herpes Zoster	Neuralgic Disorders 2	Striae after pregnancy
Ischialgie	Oedemas	Burns, 1st degree
Keloid	Pain after Preparations (grinding trauma)	Burns 1st degree, large
Lumbago chronic	Pain Prophylaxis	Manicure 1
Lumbago acute	Parodontosis	Manicure 2
Lupus Ery-thematodes	Periodontitis	Manicure 3
Luxation	Pulpa Capping	Manicure 4
Muscular Luxation	Pulpitis	Pedicure 1
Otitis Externa	Sinusitis	Pedicure 2
Paresthesia	Stomatitis	Pedicure 3
Paronchya	Super-Sensitivity of Dental Roots	Aphthae
Polyarthritits	Temporo-mandibular Joint Complaints	Eczema
Post-OP-Complaints	Trigeminus Neuralgia	Hair Loss (Alopecia Areata)
Pruritus of the Vulva	Wound Healing	Herpes Labialis
Psoriasis	WSR (Spongiosa-Therapy)	Psoriasis
Rigidity of the neck	WSR (Wound-Therapy)	Tense Muscles
Sciatica Syndrom		
Shoulder-Arm-Syndrome	applications: laser quit	applications: laser fit
Sonitus	program A	FIT craving
Tonsillitis	program B	FIT frustration
Ulcers		FIT aggression
Ulcus Vaicose		FIT eating behaviour
Vaginitis		
Wound Healing		